


TECNIS Eyhance IOL



TECNIS Eyhance IOL

Model: ICB00

DESCRIPTION	
OPTIC CHARACTERISTICS	
Powers:	+5.0 D to +34.0 D in 0.5 diopter increments
Diameter:	6.0 mm
Shape:	Biconvex, continuous, higher-order aspheric anterior surface
Material:	UV-blocking hydrophobic acrylic
Refractive Index:	1.47 at 35° C
Edge Design:	ProTEC frosted, continuous 360° posterior square edge
OPTICAL BIOMETRY*	
A-Constant (SRK/T):	119.3
AC Depth (HofferQ):	5.7 mm
Surgeon Factor (Holl.) ¹ :	1.96 mm
APPLANATION ULTRASOUND BIOMETRY	
A-Constant [†] :	118.8
Theoretical AC Depth:	5.4 mm
Surgeon Factor ¹ :	1.68 mm
HAPTIC CHARACTERISTICS	
Overall Diameter:	13.0 mm
Style:	C
Material:	UV-blocking hydrophobic acrylic
Design:	Tri-Fix , Haptics offset from optics, 1-piece lens
RECOMMENDED INSERTION INSTRUMENTS	MODEL
The UNFOLDER [®] Platinum 1 Series Screw-Style Inserter The UNFOLDER [®] Platinum 1 Series Cartridge	DK7796 1MTEC30

*Derived from clinical evaluation results of the 1-Piece IOL Platform for optical biometry.

[†]A-Constant theoretically derived for ultrasound biometry.

Reference: 1. Holladay JT. International Intraocular Lens & Implant Registry 2003. *J Cataract Refract Surg.* 2003; 29:176-197. REF2016CT0151.

The TECNIS Eyhance IOL, model ICB00, is indicated for the visual correction of aphakia in adult patients in whom a cataractous lens has been removed by extracapsular cataract extraction. The lens is intended to extend depth of focus, which improves vision for intermediate tasks, and provide similar distance vision as compared to a standard monofocal IOL. The lens is indicated for placement in the capsular bag only. For a complete listing of precautions, warnings, and adverse events, refer to the package insert.

For healthcare professionals only.

Please read the Directions for Use for Important Safety Information and consult our specialists if you have any questions.

TECNIS, TECNIS Eyhance, ProTEC, Tri-Fix and UNFOLDER are trademarks of Johnson & Johnson Surgical Vision, Inc.
©Johnson & Johnson Surgical Vision, Inc. 2018

PP2018CT5430